

Esophageal Stethoscopes

The Truer Medical disposable Esophageal Stethoscope with Temperature Sensor is used in a surgical setting for monitoring temperature as well as heart and lung sounds of anesthetized patients. Inserted in the esophagus, the stethoscope provides an accurate measurement of core body temperature as well as transmits clear heart and lung sounds to the anesthesiologist.

- Cost-effective
- Soft, thin-walled cuff provides maximum sound transmission
- Male luer allows convenient attachment to standard acoustical earpieces
- Available in 9Fr, 12Fr, 18Fr, and 24Fr to meet all size requirements
- Disposable, single-use design reduces the risk of cross-infection and ensures a quality product each time a new unit is used
- 400 series temperature sensor
- Available in both Standard and Hi-Cuff configuration
- Tube specially designed "frosted" finish to aid in the ease of insertion and removal
- Private labeling available

Ordering Information and Technical Data

Esophageal Stethoscope w/ 400 series sensor


Item No.	OD Size	Cuff Type	Box Qty	Case Qty
ES 9400	9 Fr	Standard	25	50
ES 12400	12 Fr	Standard	25	50
ES 18400	18 Fr	Standard	25	50
ES 24400	24 Fr	Standard	25	50
ES 12400H	12 Fr	Hi-Cuff	25	50
ES 18400H	18 Fr	Hi-Cuff	25	50
ES 24400H	24 Fr	Hi-Cuff	25	50

Esophageal Stethoscope w/ no sensor

ES 9000	9 Fr	Standard	25	50
ES 12000	12 Fr	Standard	25	50
ES 18000	18 Fr	Standard	25	50
ES 24000	24 Fr	Standard	25	50

Temperature Interface cables available for most monitors.

Note to Patients: All devices shown here are prescription products and must be obtained through a licensed physician or homecare provider. Patients cannot purchase directly from Truer Medical, Inc.



Hi-Cuff

Standard Cuff



400 series connector

Manufactured by Truer Medical, Inc.

1050 North Batavia, Unit C
 Orange, CA 92867, USA
 Phone: (714) 628-9785 Fax: (714) 633-3822
 www.truermedical.com

