

# General Purpose Probes

Truer Medical General Purpose Probes are designed for continuous temperature monitoring. Probes can be introduced through oral or rectal openings. Probe tips have special low-friction "frosted" finish to aid in the ease of placement and patient comfort.


- Use for routine monitoring of core body temperature
- Special "frosted" surface treatment for ease of placement and patient comfort
- Printed depth marks ensure accurate placement
- Atraumatic RF-finished tip for patient comfort
- Available in 9Fr and 14Fr sizes
- Versatile design of the general purpose probe will reduce inventory (allows for oral or rectal introduction)
- Disposable, single-use units reduce the risk of cross-infection and ensure a quality product each time a new unit is used
- 400 series temperature sensor
- Private labeling available

## Ordering Information and Technical Data

### General Purpose Probes

Item No.	Description	OD Size	Box Qty	Case Qty
GP 9400	400 series sensor	9 Fr	50	100
GP 14400	400 series sensor	14 Fr	50	100

Note to Patients: All devices shown here are prescription products and must be obtained through a licensed physician or homecare provider. Patients cannot purchase directly from Truer Medical, Inc.



GP 14400



400 series connector