

# Temperature Monitor Interface Cables

Truer Medical Temperature Monitoring Interface Cables offer several different interface configurations. By standardizing temperature probes and using an interface cable, inventory can be reduced while maintaining compatibility with most monitoring equipment.

- Cost-effective temperature monitoring solution
- Standard cables are 12' in length
- Cables are reusable
- Each cable allows connection of standard probes to different monitors
- Available with an HP, phonojack, or other standard type interfaces.
- Custom cable configurations available

## Ordering Information and Technical Data

### Temperature Monitor Interface Cables

| Item No. | Description                                     |
|----------|---|
| 70320    | 400 series with HP type connector               |
| 70321    | 400 series with right angle phonojack connector |

Note to Patients: All devices shown here are prescription products and must be obtained through a licensed physician or homecare provider. Patients cannot purchase directly from Truer Medical, Inc.



70320 - HP type connector



70321 - Right angle phonojack

## Manufactured by Truer Medical, Inc.

1050 North Batavia, Unit C  
Orange, CA 92867, USA  
Phone: (714) 628-9785 Fax: (714) 633-3822  
[www.truermedical.com](http://www.truermedical.com)

