

Tympanic Probes

The Truer Medical Tympanic Probe is designed to be a less intrusive way to accurately measure body temperature. By measuring body temperature from the tympanic membrane—through the ear canal—the patient will experience maximum comfort.

- Convenient, less-invasive alternative to esophageal or rectal probes
- Soft memory-cell foam conforms to the ear canal, reducing risk of injury to the tympanic membrane while sealing out any ambient air which could affect the temperature reading
- Disposable single-use design reduces the risk of cross-infection and ensures a quality product each time a new unit is used
- 400 series temperature sensor
- All sensor wires come coiled and taped with Co-Adhesive tape*
- Private labeling available



*Co-adhesive tape not shown

Ordering Information and Technical Data

Tympanic Probe

Item No.	Description	Case Qty
TP 400	Probe with 400 series connector	50

Also available with 48" tail.
Temperature Interface cables available for most monitors.

Note to Patients: All devices shown here are prescription products and must be obtained through a licensed physician or homecare provider. Patients cannot purchase directly from Truer Medical, Inc.



Soft memory-cell foam



400 series connector