

Gas Sampling Lines

Truer Medical Gas Sampling Lines are manufactured to the highest quality standards using medical-grade materials. Standard luer connection ensures secure fit. Bulk boxes of 500, or individually bagged in cases of 100.

- Economical
- Disposable, single patient use
- Choice of male or female luer for machine fitting
- PVC sampling lines for measurement of end-tidal CO₂
- Co-Extruded PVC/PE sampling lines available
- Available in bulk (500) or individually bagged (100/case)
- Custom sizes and configurations available
- All Gas Sampling Lines come coiled and taped with Co-Adhesive tape*
- Private labeling available
- Custom sizes available upon request



*Co-adhesive tape not shown

Ordering Information and Technical Data

STANDARD GAS LINES – In Stock

10' Co-Extruded PVC/PE

Item No.	Tube ID	Luers	Case Qty
90100	.048"	M/M	100 and 50
90101	.048"	M/F	50

10' PVC

Item No.	Tube ID	Luers	Case Qty
90150	.050"	M/M	100 and 50
90151	.050"	M/F	100



Male luer



Female luer

CUSTOM GAS LINES – Available upon request; MOQ REQUIRED.

Note to Patients: All devices shown here are prescription products and must be obtained through a licensed physician or homecare provider. Patients cannot purchase directly from Truer Medical, Inc.

Manufactured by Truer Medical, Inc.
 1050 North Batavia, Unit C
 Orange, CA 92867, USA
 Phone: (714) 628-9785 Fax: (714) 633-3822
www.truermedical.com



ISO 13485:2016